

# Romer Analytical Sampling (RAS®) Mill

Romer Labs® has a durable mill that has been specifically developed for products that are difficult to grind due to high moisture and high oil content such as pet foods, cottonseed, tree nuts or spices.



The Romer Analytical Sampling (RAS®) Mill features interchangeable grinding burrs and power drive screws to allow for the grinding of practically any product. The moving parts are made of stainless steel and flashed chrome for food safety analysis. The mill is very versatile and will easily grind samples for moisture, protein, fat, fiber, pesticide and mycotoxin analyses.

## Benefits

- *Easy-to-Clean* – Prevents sample cross-contamination
- *Powerful* – Designed for moist & high oil products
- *Safe* – Generates little heat or dust
- *Durable* – Rugged steel construction for years of use
- *Approved* – CE Approved